
# Numerical Prediction of Particles Dynamics Within a Cytometer Application to Counting and Sizing by Impedance Measurement

Damien Isèbe<sup>1</sup>

<sup>1</sup>HORIBA Medical, Parc Euromédecine, Rue du Caducée, 34184, Montpellier, FRANCE.

#### Impedance gating and particle sizing principle:

In hematology analysis, counting and sizing of blood cells are based on electrical gating in an orifice–electrode system [1]. The resulting voltage variation allows to work out particle volume measurement but various phenomenon can produce overestimations.



Ohm's law
U=R.I


microorifice

1 PARTICLE = 1 PULSE

PULSE = f(VOLUME)

Figure 1. The orifice-electrode system

Figure 2. Principle of impedance gating






Figure 3. Isovalues of the electrical field inside the orifice

Figure 4. Effect of insulation particle on the electrical field

## Numerical model for particle counting & sizing:

#### The flow model

$$\rho \frac{\partial \vec{u}}{\partial t} - \nabla \cdot \mu (\nabla \vec{u} + (\nabla \vec{u})) + \rho (\vec{u} \cdot \nabla) \vec{u} + \nabla p = 0$$

$$\nabla \cdot \vec{u} = 0$$

### The particle transport (moving boundary $\Gamma_{\rm p}$ )

$$\vec{v}\big(M_p(t),t\big) = \vec{v}_{tr}\big(M_p(t),t\big) + \vec{v}_{ang}\big(M_p(t),t\big) \times (C_p(t) - M_p(t))$$

with

$$m \cdot \frac{d\vec{v}_{tr}}{dt} = \vec{F}_{hyd}$$
 where  $\vec{F}_{hyd} = \int_{\Gamma_p} \vec{\sigma} \cdot \vec{n} \ d\Gamma_p$ 

and

$$I^* \cdot \frac{d\vec{v}_{ang}}{dt} = \vec{T} \text{ where } \begin{cases} I^* = m \cdot \int_{\Gamma_p} \left\| C_p - M_p \right\| dM_p \\ \vec{T} = \int_{\Gamma_p} (C_p - M_p) \times (\vec{\sigma} \cdot \vec{n}) dM_p \end{cases}$$

#### The electrical field

$$-\nabla \cdot \left[ \left( \sigma_e + \frac{\epsilon}{\tau} \right) \nabla V \right] = 0$$

#### The voltage pulse

$$U_p(t) = \frac{W(t_2) - W(t_1)}{(\Delta t) \cdot \overline{I}} \text{ where } t = \frac{t_1 + t_2}{2} \text{ and } W = \frac{\epsilon}{2} \int_{\Omega} \|\vec{E}\|^2$$

**Computational Methods**: Deformed mesh and FSI setting – ALE formulation for the moving mesh. Automatic remeshing (solver stops when mesh quality is lower than 0.5) – Around 30 remesh/start cycles.




Figure 5. Remesh/start cycle to control mesh quality



Figure 6. The whole displacement of the particle through the orifice

**Results**: Considered particles are shaped like capsules and supposed to be rigid and totally insulating.




Figure 8. Velocity field – Particle through the micro-orifice (FSI)

# Effect of particle size on voltage pulses

# 

Figure 9. Effect of particle size

Figure 10. Effect of particle orientation

Effect of particle orientation

### Effect of particle trajectory on voltage pulses




Figure 11. Effect of particle trajectory – (left) numerical; (right) experimental

## Role of hydrodynamic focusing in cell sizing




Figure 12. Improvement of Particle Size Distribution by Sample Hydrodynamic Focusing

**Conclusions**: A fully innovative numerical approach, in lines with experimental tests, is presented to tackle the problem of counting and sizing particles. Substantial improvements on the transport model, including 3D processing and deformability of particles under strong hydrodynamic stresses, are investigated.

#### References:

- 1. V. Kachel and al., Elecrical Resistance Pulse Sizing: Coulter Sizing, Flow Cytometry & Sorting, 35p., 1990
- 2. D. Isèbe and al., Numerical simulation of particle dynamics in an orifice-electrode systeme, IJNMBE, 13p., 2013

